
Automated Micromanipulator (AM)

Automated Micromanipulator

Instruction Manual (Translation of the original instructions)

Software Version: PreSens Profiling Studio version 1.0.0 **Document filename:** IM_AM_Translations Original_dv2

1 Safety Notes

This device is only intended for use in the laboratory by qualified personnel according to this instruction manual and the following safety guidelines:

- ! Before operating the Automated Micromanipulator (AM), please make sure that the travel path of your microsensor is free and unobstructed by any object other than your sample. Also, make sure that the programmed traveling path of the microsensor never exceeds the size of your vessel or the distance to the underlying surface or base plate.
- In case of an emergency pull the red power plug for immediate emergency stop.
- ! The potentiometric control knob for manual movement has to be in zero position before putting the AM back into operation after an emergency stop.
- ! Always use the delivered protective bag or operate the AM in a safe environment (such as a fume hood) to avoid any harm by possible breakage of the microsensor (e. g. in samples with unknown or inhomogeneous composition or consistency).
- Never operate the manual movement of the linear stage without direct supervision.
- ! The device may produce enough force to cause personal injury. Be careful to keep hair, body parts, jewellery, and clothing from being caught in moving components. Warning labels are used on the devices to indicate areas of particular concern.
- During continuous operation, the device's motor may feel hot to the touch. Although this is normal, care should be taken when handling the device. If the device emits a burnt smell, it may be damaged, in which case you should immediately cut the power supply, cease operation and

2 Description of the Automated Micromanipulator

Fig. 1 Automated Micromanipulator mounted to the Heavy Stand

The Automated Micromanipulator (AM) is specifically designed for profiling applications with the PreSens profiling microsensors (PM), and can also be operated with needle-type housed (NTH) and implantable (IMP) microsensors. The system allows moving the microsensor vibration-free with μ m reading accuracy and enables exact localization of the sensor in the sample. Automated profiling can be performed along one dimension over a maximum travel range of 7.5 cm in step sizes down to 10 μ m. The micromanipulator additionally comprises a tilting platform, which allows adjusting the microsensor at an angled position. The associated, user-friendly, and database-supported software PreSens Profiling Studio enables control of the Automated Micromanipulator and the respective oxygen or pH meter. The AM is compatible with PreSens oxygen meters Microx TX3 (sensor type PSt1), Microx TX3 trace (with PSt1 sensor type only), Microx 4 (sensor type PSt7), Microx 4 trace (sensor type PSt7 & PSt8) and the pH meter pH-1 micro (sensor type HP5). The software offers multiple features from clear data organization and export, easy creation of profiling templates, to initial analysis functions.

Features

- Fully automated system
- Software 'PreSens Profiling Studio' included
- Plug & play USB-connection
- Manual Micromanipulation in 3D
- O Automated profiling with 1 μm resolution
- For all PreSens microsensors and micro transmitters

3 Use of the Automated Micromanipulator

- 1. Check your shipment of the AM for completeness:
 - 1 x Manual micromanipulator
 - 1 x Linear stage with attached tilting platform
 - 1 x Profiling Studio software (CD)
 - 1 x Power adapter
 - 0 1 x USB to serial-mini cable
 - 1 x 7-port USB 2.0 hub
 - 3 x Country adapters (UK, US, AU)
 - Install the PreSens Profiling Studio software on your PC / notebook. Follow the onscreen instructions of the installer.
 - Fasten the micromanipulator to the linear stage with the two Allen-socket screws using the Allen key.
 - 4. Loosen the locking screws of the tilting platform, tilt the AM and place it sideways so the knurled screw fits through the drill hole underneath. Attach the screw to a nut block or M6 thread in your micromanipulator stand and tighten it with the respective key.
 - 5. Insert the red power adapter plug in the power connector on the AM and connect the power adapter to the power grid. Connect the USB cable of the AM to the USB hub delivered with the AM. Connect the USB hub to your PC / notebook and its power adapter to the power grid.
 - Connect the microsensor to the sensor connector on the oxygen / pH meter and mount it in the sensor holder of the AM.
 - Connect the transmitter to the power supply (not necessary with Microx 4 / Microx 4 trace) and to a USB port of the USB hub connected to the PC / notebook. Switch on the transmitter.
 - Close all other applications as they might interfere with the software and start the PreSens Profiling Studio.

Power Supply

Meter

- 1 x Power strip with overvoltage protection
- 1 x Clamp for bare fiber microsensors
- 1 x Knurled screw (M6)
- 2 x Allen-socket screws (M6)
- 1 x Allen key (size 4.0)
- 1 x Key for knurled screw

Fig. 2: Attach the micromanipulator to the linear stage.

Fig. 3: Mount the AM to your micromanipulator stand.

Fig. 4: Connect the AM to the power supply (left) and a PC / notebook (right).


Fig. 5: Connect the microsensor and mount it on the AM.

4 Technical Data

Specifications

compatibility	Profiling (PM), Needle-type housed (NTH) and Implantable (IMP) oxygen & pH microsensors
ravel range automated	x-axis: 75 mm
ravel range manual	x-axis: 37 mm, fine drive 10 mm
J	y-axis: 20 mm
	z-axis: 25 mm
lax. motor velocity	4 mm/s
lin. motor velocity	1 μm/s
esolution	1 μm
epeatability	< 2.5 µm
laximum continuous thrust	100 N
laximum cantilever load	100 N cm
lounting adapter	M6 screw, 13 mm length
ower Supply	
upply voltage	15 VDC (2.1 mm center positive plug)
ower adapter	100 – 240 VAC, 50 – 60 Hz (Use supplied power adapter only)
nvironmental Conditions	
perating temperature	0 to +50 °C
lax. relative humidity	Up to 80 % (non-condensing)
imensions / Weight	
imensions	275 mm x 95 mm x 220 mm
Veight	2070 q
<u> </u>	•
igital Interface	
ISB interface cable to PC	Cable included

5 Operational Notes

Storage

Before storing your device, retract any extended components in order to keep them clean and to protect them from damage. Store the device in a dry environment.

Environmental Conditions

The device is NOT splash or water proof. Do not submerse the device in liquid. Do not expose the device to vibration or shock, or extreme conditions, such as temperatures exceeding device ratings, radiation, and dusty or humid environments.

Power Adapter

The Automated Micromanipulator always has to be used with the original power adapter (Input: 100/240 VAC, 50-60 Hz; Output: 15 VDC / 1.2 A, max. 16.2 W, type 3A-183WP15) which is supplied.

Maintenance

The Automated Micromanipulator is maintenance-free.

The housing should be cleaned with a cloth only. Avoid any moisture

entering the housing. Never use benzine, acetone, alcohol or any other organic solvents for cleaning.

Service

Alignment, rework or repair work may only be carried out by the manufacturer:

PreSens

Precision Sensing GmbH Am Biopark 11

93053 Regensburg Germany

Phone +49 941 94272100 Fax +49 941 94272111

info@PreSens.de

Please contact our service team in case of any questions. We look forward to helping you and are open for any proposition or criticism.

6 CE Conformity

The equipment is confirmed to comply with the requirements set out in the Council Directive relating to Electromagnetic Compatibility (2004/108/EEC), Low Voltage (2006/95/EEC) and machinery (2006/42/EC). For the evaluation of the above mentioned Council Directives the following harmonized standards were consulted:

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PreSens Precision Sensing GmbH