

Vial Adapter VA-20mL-Single



Do not use excessive force when tightening the screws of the Vial Adapter! Applying too much force might break the glass vials!



Do not use other screws than the ones supplied! Too large screws might break the SensorVial and cause severe harm.

Please contact our service team at info@presens.de in case you need new screws.



Do not use other glass vials with the Vial Adapter! Too large glass vials might break and cause severe harm.

Please contact our service team at info@presens.de in case your application requires different vial formats.

The Vial Adapter VA-20mL-Single is for use with the 20 mL-SensorVials with integrated sensor foil. It consists of a plastic ring with 3 holes, and 3 slotted head plastic screws. The 3 screws are identical.



Fig. 2 Adapter ring



Fig. 1 Adapter screws

Two screws are placed at the sides for attaching the Vial Adapter to the SensorVials. The screw from the top holds a light guide for read-out of the sensor in place. Fix the screws with a screwdriver.



Fig. 3 Vial Adapter with attached screws



Fig. 4 Vial Adapter attached to SensorVial